


General Purpose, Simple to Use Polarimeter

sugar
pharmaceutical
chemical
education
cosmetics
flavours
fragrances
essences



Bellingham + Stanley Inc.
Quality Scientific Instruments

ADP410 Polarimeter

Featuring:

- Dual scale
- Temperature display
- Optical Density display
- Temperature compensation
- LED/Interference filter
- RS232 (print/PC)
- Bright LED display
- Easy operation

Optional Accessories:

- Flow tubes
- Water-jacketed tubes
- Slotted lids
- Quartz control plates
- Printers

The ADP410 is a dual scale, fully automatic polarimeter designed for use in many applications that require measurement of optical rotation.

The instrument is housed in a rugged chemical resistant case making it suitable for use in factory environments as well as in the laboratory. Standard, jacketed and flow type tubes may be used, possibly requiring the use of specially suited slotted lids.

Operation is very simple by way of four graphically identified push buttons including scale and compensation selection, zero calibration and the displaying of optical density and temperature. Results are shown on a bright LED display and may be printed via an RS232 interface.

Temperature compensation can be achieved from a single sensor, which can measure either the temperature within the sample chamber or, when placed in the center-filling arm of the measurement tube, the temperature of the actual sample.

The sample reading value and optical density are alternately displayed by repeatedly pressing the READ key. The instrument can be checked and calibrated using quartz control plates by accessing the span facility in the setup menu. If a sample is too dark, the instrument will display an error message.

An LED and narrow band interference filter provides illumination of the sample and requires no maintenance.



ADP410 instrument with Polarimeter tube in place.



ADP410 instrument display panel.

Specifications

	Angular (°A)	ISS (°Z)
Range	-90 to +90	-225 to +225
Resolution	0.01	0.01
Accuracy	±0.02	±0.05

Light Source	LED/Interference filter (589.3nm)
Beam Diameter	4 mm
Max. Tube length	220 mm
Reading Time	15 seconds
Temperature	5-40 °C
Temperature Compensation	None, sugar, quartz
Optical Density Range	0.0 to 2.0 OD

Order Code Description

37-10	ADP410 digital polarimeter supplied with one center filling tube (200mm), instruction manual and certificate of conformity.
34-20	Quartz control plate 34 °A / 100 °Z
34-21	Quartz control plate 5 °A / 15 °Z
37-009	Slotted lid for 200mm center fill tube with water jacket
37-010	Slotted lid for 100mm center fill tube with water jacket
37-011	Slotted lid for 200mm flow type tube with water jacket
37-012	Slotted lid for 100mm flow type tube with water jacket

Nett Weight: 9.0 kg

Actual Dimension:

620(W) x 340(D) x 150(H) mm

Power Supply: External 90-250v~, 50/60 Hz.
Supplied with instrument.

Bellingham+Stanley Ltd pursues a policy of continuous product development and improvement, and, as such, information given on this Data Sheet may be updated or withdrawn without notice.



Bellingham + Stanley Inc.

1000 Hurricane Shoals Road, Building D, Suite 300, Lawrenceville, GA 30043 USA

Phone: 770 822 6898. Fax: 770 822 9165

www.bellinghamandstanley.com sales@bs-rfm-inc.com

1201/0305