


# Printers & Interface Cables

The direct recording of data from scientific instruments is an important requirement in quality control. A range of printers is available for use with Bellingham + Stanley refractometers and polarimeters, depending on the application and model of instrument. For users adopting Good Laboratory Practice and with instruments that have a real-time clock, the CBM-910 is the recommended printer.

For users requiring a GLP record from an instrument that does not support date/time printing, the DP-24 may be used, as it incorporates its own clock.



CBM-910  
Miniature Dot  
Matrix Printer



DP-24 Desk  
Top Printer



## CBM-910 Miniature Dot Matrix Printer

- Miniature 24-column dot matrix impact printer.
- Serial (RS232C) or Parallel (Centronics) models available.
- Ideal for RFM300+, RFM700 and RFM900 digital refractometers and all ADP polarimeters.
- Recommended for Good Laboratory Practice applications.
- Supplied with Euro/UK or US mains adapter as indicated.

## DP-24 Desk Top Printer\*

- Miniature 24-column dot matrix printer.
- Integral date/time stamp.
- Recommended for GLP applications where clock is not available in the instrument.
- 90-230V, 50-60 Hz mains.
- Supplied with mains cable.

\* Not suitable for damp environments unless suitably protected

Order Code	Description	Model
55-11 55-13	CBM-910 Printer Parallel: Euro/UK - 230V <b>CBM-910 Printer Parallel: US - 110V</b>	RFM300+/RFM900 & ADP440
55-14 55-16	CBM-910 Printer Serial: Euro/UK - 230V <b>CBM-910 Printer Serial: US - 110V</b>	All RFM and ADP instruments
55-17	DP-24 Printer Serial (GLP: Date/Time) 90 to 230V	RFM700/ADP410 & ADS420

NOTE: Interface cables are not supplied. Please select from below:

## Accessories

Order Code	Description	Model
54-01	Printer Cable for CBM-910 Serial	RFM100/RFM300
54-02	Printer Cable for CBM-910 or BS1012 Serial	RFM900/RFM700 & ADP/ADS
54-03	Printer Cable for CBM-910 Parallel	RFM300+/RFM900/ADP440
54-08	Printer Cable for BS1012 Parallel	
54-09	Printer Cable for DP-24	RFM700/ADP410 & ADS420
55-91	Paper Roll (Qty 20)	CBM-910/DP-24 & BS1012

## Accessories - Instrument to PC Interface Cables

Order Code	Description
54-01 & 54-51 (set)	RFM100 & RFM300
54-07	RFM700/900 & all ADP/ADS

### Head Office

Longfield Road, Tunbridge Wells,  
Kent TN2 3EY, United Kingdom  
Phone: +44 (0) 1892 500400  
Fax: +44 (0) 1892 543115  
sales@bellinghamandstanley.co.uk

### USA Office

1000 Hurricane Shoals Road, Building D,  
Suite 300, Lawrenceville, GA 30043 USA  
Phone: 770 822 6898  
Fax: 770 822 9165  
sales@bs-rfm-inc.com

[www.bellinghamandstanley.com](http://www.bellinghamandstanley.com)



**Bellingham  
+ Stanley**